

# Plug in Gait WebEx Training Session 3

## *Interpreting PiG results: PiG biomechanical modelling*

Gabriele Paolini  
Support Engineer

## *INTRODUCTION*

What is Plug in Gait??

## INTRODUCTION

### What is Plug in Gait??

Plug in Gait is the commercial name used by Vicon for the implementation of what is commonly called 'Conventional Gait Model' (CGM)

The CGM is a biomechanical model for the lower limbs developed by Kadaba, Davis and the Helen Hayes Hospital.


Kadaba et al., 1989, Journal of Orthopaedic Research.

Kadaba et al., 1990, Journal of Orthopaedic Research.


Davis et al., 1991, Human Movement Sciences.

## WHAT IS A MODEL?


IN GENERAL, A MODEL IS A MATHEMATICAL RELATION THAT ALLOWS FORESEEING THE BEHAVIOR OF THE SYSTEM THAT HAS BEEN MODELLED, GIVEN SOME KNOWN INPUTS AND PARAMETERS



*BIOMECHANICAL MODELLING AND PLUG IN GAIT*



## BIOMECHANICAL MODELLING AND PLUG IN GAIT



Plug In Gait contains four modelling modules:

- |                               |   |                                 |
|-------------------------------|---|---------------------------------|
| 1. Lower Body kinematic model | → | Rigid body segments definition  |
| 2. Upper Body kinematic model |   | Joint Angles, Joint Centers     |
| 3. Lower Body kinetic model   |   | Mass and Inertia Integration    |
| 4. Upper Body kinetic model   | → | Joint Forces, Moments and Power |

## *PLUG IN GAIT GENERAL KINEMATIC MODELLING CONSIDERATIONS*

1. Rigid Body Hypothesis
2. Rigid segments are defined starting from
  1. Physical markers
  2. Virtual markers calculated using physical markers and subject measurements
3. Rigid segments are defined on a frame-to-frame basis



It is vitally important to place the markers accurately and to get good and stable 3D marker positions reconstructions

## *PLUG IN GAIT GENERAL KINETIC MODELLING CONSIDERATIONS*

1. Kinetic calculations are always performed alongside kinematics
2. If Force Platforms (FPs) are present, the Ground Reaction Forces (GRFs) are added automatically to the processing workflow
3. Kinetic modelling is performed even when no FPs are present
  1. In this case, the kinetic calculations are only valid during swing phase
4. The anthropometric measures – segment masses, centres of gravity and radii of gyration – are taken from David Winter's tables (add ref)
5. Joint kinetics are calculated using the inverse dynamics procedure




## *PLUG IN GAIT SUBJECT MEASUREMENTS*

Subject Measurements represent the parameters of the model

## PLUG IN GAIT SUBJECT MEASUREMENTS

Subject Measurements represent the parameters of the model



## *PLUG IN GAIT SUBJECT MEASUREMENTS*

Subject Measurements represent the parameters of the model

### REQUIRED SUBJECT MEASUREMENTS

#### General

*Body Mass*

*Height*

#### Lower Body

*Leg Length*

*Knee Width*

*Ankle Width*

#### Upper Limbs

*Shoulder Offset*

*Elbow Width*

*Wrist Width*

*Hand Thickness*

## *PLUG IN GAIT INPUT*

1. A minimum number of markers need to be present in the trial
  1. Pelvis markers for lower body model
  2. Thorax markers for upper body model
2. Alternative marker sets possible
  1. SACR marker or PSI markers
  2. Wand markers for SACR, THI, TIB
  3. KADs
  4. Medial malleolus

EACH OF THE ABOVE ALTERNATIVES IS VALID AND NO SPECIAL ACTION NEEDS TO BE TAKEN TO RUN IT.

## PLUG IN GAIT OUTPUT

The outputs generated by PiG are saved as trajectories in the c3d file

Each PiG output is in fact a vectorial physical entity – i.e. is composed by three components X, Y, Z – the same as the trajectory of a physical marker

A message regarding the successful or unsuccessful running of PiG appears in the Nexus processing log, by clicking on the *Log* tab in the *Communications* window

PiG outputs are visible within Nexus, under the subject data tree in the *Resources* tools pane

## KINEMATICS

### 1. Joint Angles

- a. Relative angles between two rigid segments, i.e. joint angles. Always ordered as Flexion, Abduction, Rotation.
- b. Absolute angles between a rigid segment and the laboratory fixed reference. Always ordered as rotations around the global X, Y and Z axes

### 2. Bones

- a. A set of four virtual points will be associated to each modelled body segment
- b. *segName#O*: Origin of the *segName* segment
- c. *segName#L*: Direction of the lateral axis of the *segName* segment
- d. *segName#A*: Direction of the anterior axis of the *segName* segment
- e. *segName#P*: Direction of the vertical axis of the *segName* segment
- f. The Bones output section is mainly used for positioning the bone meshes in Polygon

# KINETICS

## 1. Joint Forces

- a. Net Joint Forces expressed in the local reference system of each rigid body segment
- b. Units: [N/Kg]

## 2. Joint Moments

- a. Net Joint Moments estimated solving the equations of motion for the segments of the lower limbs – excluding the Pelvis (Ramakrishnan et al, 1987)
- b. Anthropometric Measurements from published tables (Winter)
- c. 'External Forces' convention is used: a GRF that would cause extension produces a positive extension moment
- d. Net Joint Moments ordered as Flexion, Abduction, Rotation
- e. Units: [Nmm/Kg]

## 3. Joint Powers

1. Scalar product between Joint Moments and Joint Angular Velocities
2. Joint Powers can be expressed as a scalar or as three components
3. Units: [W/Kg]

***PLUG IN GAIT LOWER BODY  
KINEMATIC MODELLING***



## *PLUG IN GAIT LOWER BODY KINEMATIC MODELLING*

Let's revise the main hypotheses behind PiG kinematic modelling:

1. Rigid Body Hypothesis
2. Rigid segments are defined starting from
  1. Physical markers
  2. Virtual markers calculated using physical markers and subject measurements
3. Rigid segments are defined on a frame-to-frame basis




It is vitally important to place the markers accurately and to get good and stable 3D marker positions reconstructions

## *PLUG IN GAIT LOWER BODY KINEMATIC MODELLING*


The lower body PiG kinematic model includes 7 rigid segments:

1. Pelvis
2. Left Thigh
3. Right Thigh
4. Left Shank
5. Right Shank
6. Left Foot
7. Right Foot

Pelvis



## PELVIS MARKERS




NOTE: Using the two posterior *PSI* markers or the *SACR* marker is irrelevant.

When the *SACR* marker is present, it will be used for modelling.


When the *PSI* markers are present, the mid point between them will be calculated and used for modelling. If only one *PSI* marker is visible, that will be used for modelling

*PELVIS TECHNICAL REFERENCE SYSTEM*




## PELVIS TECHNICAL REFERENCE SYSTEM

- Origin:  $(RASI + LASI)/2$



## PELVIS TECHNICAL REFERENCE SYSTEM


- Origin:  $(RASI + LASI)/2$
- Y axis direction:  $RASI \rightarrow LASI$





## PELVIS TECHNICAL REFERENCE SYSTEM


- Origin:  $(RASI + LASI)/2$
- Y axis direction:  $RASI \rightarrow LASI$
- Z axis direction: Perpendicular to the plane defined by  $LASI, RASI, SACR$






## PELVIS TECHNICAL REFERENCE SYSTEM

- Origin:  $(RASI + LASI)/2$
- Y axis direction:  $RASI \rightarrow LASI$
- Z axis direction: Perpendicular to the plane defined by  $LASI, RASI, SACR$



## PELVIS TECHNICAL REFERENCE SYSTEM

- Origin:  $(RASI + LASI)/2$
- Y axis direction:  $RASI \rightarrow LASI$
- Z axis direction: Perpendicular to the plane defined by  $LASI, RASI, SACR$
- X axis direction: cross product between Y and Z unit vectors



## HIP JOINT CENTRES CALCULATION

- Newington – Gage model used for defining the positions of the Hip Joint Centres (HJCs) in the pelvis technical reference system (Davis et al, 1991)
- The Newington – Gage model uses the Asis to Trochanter distance and the Inter-Asis distance to create the position vector of the HJCs in the Pelvis technical reference system
- Asis to Trochanter and Inter-Asis distances can be manually added to the subject measurements (optional measurements)
- Asis to Trochanter and Inter-Asis distances can be automatically calculated by PiG, if not present in the subject measurements
  - Inter-Asis = DIST(LASI,RASI)
  - AsisTrocDistance =  $0.1288 * \text{LegLength} - 48.56$

## HIP JOINT CENTRES CALCULATION

$$LHJC_x = C * \cos(\mathcal{G}) * \sin(\beta) - (AsisTrocDist + mm) * \cos(\beta)$$

$$LHJC_y = -(C * \sin(\mathcal{G}) - aa)$$

$$LHJC_z = -C * \cos(\mathcal{G}) * \cos(\beta) - (AsisTrocDist + mm) * \sin(\beta)$$

Where:

$$\mathcal{G} = 0.5rad$$


$$\beta = 0.314rad$$

$$AsisTrocDist = 0.1288 * LegLength - 48.56 *$$

$$C = MeanLegLength * 0.115 - 15.3$$

$$aa = (InterAsis) / 2$$

$mm$  = marker radius




\*: It can also be entered in the Subject Measurements

## PELVIS ANATOMICAL REFERENCE SYSTEM

Once the HJCs are defined, the origin of the reference system associated to the Pelvis segment is shifted to the mid point of the HJCs.


The axes orientation remains unchanged



## PELVIS ANATOMICAL REFERENCE SYSTEM

Once the HJCs are defined, the origin of the reference system associated to the Pelvis segment is shifted to the mid point of the HJCs.

The axes orientation remains unchanged





## PELVIS SUMMARY

### ☺ Markers

☺ *LASI, RASI, LPSI, RPSI*

OR

☺ *LASI, RASI, SACR*

### ☺ Technical reference definition

### ☺ Hip Joint Centres calculation

☺ InterAxis distance from subject measurements OR calculated by PiG

☺ ASIS to Trochanter distance from subject measurements OR calculated by PiG

### ☺ Anatomical reference definition

☺ Origin:  $(RHCJ + LHJC)/2$

☺ Y axis direction:  $RASI \rightarrow LASI$

☺ Z axis direction: Perpendicular to the plane defined by *LASI, RASI, SACR*

☺ X axis direction: cross product between Y and Z unit vectors

Thigh






## THIGH MARKERS

- **KNE:** On the most lateral aspect of the femur lateral condyle
- **THI:** On the lateral aspect of the thigh lying on the femoral frontal plane – defined by *HJC*, lateral and medial femoral condyles. Height not critical.

NOTE: The *THI* marker placement is very important since, together with the *KNE* marker and the Hip Joint Center, defines the orientation in space of the femoral frontal plane



## THIGH ANATOMICAL REFERENCE SYSTEM

By looking at the marker set for the thigh we notice that there are only two physical markers attached to the thigh segment.


To describe a rigid body motion in space at least three points are needed.

*HJC?*

Yes, but not only.

PiG in fact calculates the location of the Knee Joint Center (*KJC*) so that the direction of the longitudinal axis of the femur can be defined.

The *KJC* position is determined by applying what we call the CHORD function.




## KNEE JOINT CENTER CALCULATION

The *KJC* is defined as the point at distance *KneeOS* from the *KNE* marker in the plane defined by *KNE*, *THI* and *HJC*.

The angle *KNE-KJC-HJC* must be 90°


$KneeOS = (MarkerDiameter + KneeWidth) / 2$

*KneeWidth* from Subject Measurements (**Required**)




## THIGH ANATOMICAL REFERENCE SYSTEM

- Origin: *KJC*




## THIGH ANATOMICAL REFERENCE SYSTEM

- Origin: *KJC*
- Z axis direction: *KJC* → *HJC*




## THIGH ANATOMICAL REFERENCE SYSTEM

- Origin: *KJC*
- Z axis direction: *KJC* → *HJC*
- X axis direction: perpendicular to the plane defined by *HJC*, *KNE*, *THI*



## THIGH ANATOMICAL REFERENCE SYSTEM


- Origin: *KJC*
- Z axis direction: *KJC* → *HJC*
- X axis direction: perpendicular to the plane defined by *HJC*, *KNE*, *THI*





## THIGH ANATOMICAL REFERENCE SYSTEM

- Origin: *KJC*
- Z axis direction: *KJC*  $\rightarrow$  *HJC*
- X axis direction: perpendicular to the plane defined by *HJC*, *KNE*, *THI*
- Y axis direction: cross product between Z and X unit vectors





## THIGH SUMMARY

### ☺ Markers

☺ *KNE, THI*

☺ The anterior-posterior position of the *THI* marker sets the orientation of the femoral frontal plane

### ☺ Knee Joint Center calculation

☺ The CHORD function creates a virtual point (*KJC*) at a distance *KneeOS* from the *KNE* marker, on the plane defined by *KNE, THI, HJC*

☺ The *KNE-KJC-HJC* angle is  $90^\circ$

### ☺ Anatomical reference definition

☺ Origin: *KJC*

☺ Z axis direction: *KJC*  $\rightarrow$  *HJC* (up)

☺ X axis direction: Perpendicular to the plane defined by *HJC, KNE, THI* (anterior)

☺ Y axis direction: cross product between Z and X unit vectors (towards the left for both left and right limbs)


Shank



## SHANK MARKERS

- **TIB**: On the lateral aspect of the shank lying on the tibial frontal plane – defined by *KJC*, lateral and medial malleoli. Height not critical.
- **ANK**: On the most lateral aspect of the lateral malleolus

NOTE: The *TIB* marker placement is very important since, together with the *KJC* and the *ANK* marker, defines the orientation in space of the tibial frontal plane



## SHANK ANATOMICAL REFERENCE SYSTEM

By looking at the marker set for the shank we notice that there are only two physical markers attached to the segment.


To describe a rigid body motion in space at least three points are needed.

*KJC?*

Yes, but not only.

PiG in fact calculates the location of the Ankle Joint Center (*AJC*) so that the direction of the longitudinal axis of the femur can be defined.

The *AJC* position is determined by applying the CHORD function, similarly to the thigh segment




## ANKLE JOINT CENTER CALCULATION

The *AJC* is defined as the point at distance *AnkleOS* from the *ANK* marker in the plane defined by *ANK*, *TIB* and *KJC*.

The angle *ANK-AJC-KJC* must be 90°

$AnkleOS = (MarkerDiameter + AnkleWidth)/2$

*AnkleWidth* from Subject Measurements  
**(Required)**



## SHANK ANATOMICAL REFERENCES

PiG creates two anatomical shank references.

### • TORSIONED TIBIA

- Tibial frontal plane being defined by *ANK*, *TIB* and *KJC*
- Used to describe the distal motion of the shank complex (*ankle angles*)


### • UNTORSIONED TIBIA

- Torsioned Tibia segment rotated around its z axis by *Tibial Torsion* degrees
- The *Tibial Torsion* is entered in the subject measurements (**optional**)
- The direction of the segment's frontal plane is parallel to the knee flexion axis
- Used to describe the proximal motion of the shank complex (*knee angles*)

## SHANK ANATOMICAL REFERENCE SYSTEM

### TORSIONED TIBIA


- Origin: AJC



## SHANK ANATOMICAL REFERENCE SYSTEM

### TORSIONED TIBIA

- Origin: AJC
- Z Axis direction: AJC → KJC






## SHANK ANATOMICAL REFERENCE SYSTEM

### TORSIONED TIBIA


- Origin: AJC
- Z Axis direction: AJC → KJC
- X axis direction: perpendicular to the plane formed by TIB, AJC, KJC



## SHANK ANATOMICAL REFERENCE SYSTEM

### TORSIONED TIBIA


- Origin: AJC
- Z Axis direction: AJC → KJC
- X axis direction: perpendicular to the plane formed by TIB, AJC, KJC



## SHANK ANATOMICAL REFERENCE SYSTEM

### TORSIONED TIBIA

- Origin: AJC
- Z Axis direction: AJC  $\rightarrow$  KJC
- X axis direction: perpendicular to the plane formed by TIB, AJC, KJC
- Y axis direction: cross product between Z and X unit vectors




## SHANK ANATOMICAL REFERENCE SYSTEM

### UNTORSIONED TIBIA

- Origin: AJC
- Z Axis direction: AJC → KJC
- X axis direction: perpendicular to the plane formed by THI, KJC, HJC – **same as Thigh segment!!!**
- Y axis direction: cross product between Z and X unit vectors

NOTE: To obtain the untorsioned tibia, Plug In Gait rotates the torsioned tibia around its Z axis by *Tibial Torsion* degrees (from the subject measurements)



## SHANK SUMMARY

### ☺ Markers

- ☺ *TIB, ANK*

### ☺ Ankle Joint Center calculation


- ☺ *AJC* at a distance *AnkleOS* from the *ANK* marker, on the plane defined by *ANK, TIB, KJC*. The angle *ANK-AJC-KJC* measures  $90^\circ$

### ☺ Shank Segments Definition

- ☺ Torsioned Tibia
- ☺ Untorsioned Tibia

### ☺ Anatomical reference definition

- ☺ Origin: *AJC*
- ☺ Z axis direction: *AJC* → *KJC* (up)
- ☺ X axis direction: Perpendicular to the plane defined by *ANK, TIB, AJC* (anterior)
- ☺ Y axis direction: cross product between Z and X unit vectors (left)




Foot


## FOOT MARKERS

- **TOE:** On the second metatarsal head, on the mid-foot side of the equinus break between forefoot and mid-foot
- **HEE:** On the calcaneus. The imaginary line connecting HEE and TOE should be parallel to the longitudinal axis of the foot

NOTE: The *TOE* and *HEE* markers should be at the same height from the foot plantar surface. However, for consistency, this can be achieved by setting the 'assume foot flat' option in PiG




# FOOT ANATOMICAL REFERENCE SYSTEM - WORKFLOW






## FOOT ANATOMICAL REFERENCE SYSTEM - WORKFLOW




## FOOT UNCORRECTED ANATOMICAL REFERENCE SYSTEM

- Origin: *TOE*
- Z Axis direction: *TOE* → *AJC*




## FOOT UNCORRECTED ANATOMICAL REFERENCE SYSTEM

- Origin: *TOE*
- Z Axis direction: *TOE* → *AJC*
- Y axis direction: perpendicular to the plane formed by *TOE*, *AJC*, *KJC*




## FOOT UNCORRECTED ANATOMICAL REFERENCE SYSTEM

- Origin: *TOE*
- Z Axis direction: *TOE* → *AJC*
- Y axis direction: perpendicular to the plane formed by *TOE*, *AJC*, *KJC*



## FOOT UNCORRECTED ANATOMICAL REFERENCE SYSTEM

- Origin: *TOE*
- Z Axis direction: *TOE* → *AJC*
- Y axis direction: perpendicular to the plane formed by *TOE*, *AJC*, *KJC*
- X axis direction: cross product between Y and Z unit vectors





FOOT STATIC REFERENCE

YES

Foot flat option?

NO





FOOT STATIC REFERENCE

YES

Foot flat option?

NO





# FOOT STATIC REFERENCE

YES

Foot flat option?

NO






For dynamic trials, the static offsets are applied to the foot anatomical reference:

$\alpha$  : Static Plantar-Flexion Offset

$\beta$  : Static Rotation Offset



## FOOT SUMMARY

### ☺ Markers

☺ *HEE, TOE*

### ☺ Anatomical reference definition

☺ Origin: *TOE*

☺ Z Axis direction: *TOE* → *AJC*

☺ Y axis direction: perpendicular to the plane formed by *TOE, AJC, KJC*

☺ X axis direction: cross product between Y and Z unit vectors

### ☺ If the trial is static, foot static reference definition

☺ Depends on foot flat option

☺ Depends on the Sole Delta measurement

☺ Plantar-Flexion and Rotation Offset calculation

☺ If the trial is dynamic, the Plantar-Flexion and Rotation offsets are applied to the foot anatomical reference


## REFERENCES

Kadaba, M.P., Ramakrishnan, H.K., Wootten, M.E, Gainey, J., Gorton, G. & Cochran, G.V.B (1989). *Repeatability of kinematic, kinetics and electromyographic data in normal adult gait*. Journal of Orthopaedic Research, 7, 849-860.

Kadaba, M.P., Ramakrishnan, H.K. & Wooten, M.E. (1990). *Lower extremity kinematics during level walking*. Journal of Orthopaedic Research, 8, 849-860.

Davis, R., Ounpuu, S., Tyburski, D. & Gage, J. (1991). *A gait analysis data collection and reduction technique*. Human Movement Sciences, 10, 575-587.

# PLUG IN GAIT WORKFLOW




Subject Action  
(in Capture Volume)

Operator Action

Nexus /Plug-in Gait  
Action

5



6

